




DOUBLE MAGNETIC STIRRER

MODEL DMS • CODE 1776

QUANTITY	CONTENTS	CODE
1	Magnetic Stirrer	26502
4	Magnetic Stirring Bars	2-2185
1	AC Adapter, 9V	1744

To order replacement components, use the specified code number.



INTRODUCTION

The Magnetic Stirrer will hold two 250 mL flasks or beakers. Place containers on stirrer surface and turn stirrer on. Adjust the stirring speed with the individual control knobs.

PROCEDURE

1. Turn both speed control knobs fully counterclockwise.
2. Plug adapter into side of stirrer. Connect to appropriate power source.
3. Place flat bottomed flask containing solution on stirrer.
4. Place stirring bar in flask. Move flask so stirring bar is centered.
5. Press "ON/OFF" button to turn stirrer on. Red indicator light should be on.
6. Adjust speed using speed control knob. Turn knob clockwise to increase stirring; turn knob counterclockwise to slow stirring.
7. After mixing is complete, turn speed control knob fully counterclockwise. Push "ON/OFF" button to turn stirrer off. Red indicator light should be off.

LaMOTTE COMPANY

Helping People Solve Analytical Challenges®

PO Box 329 • Chestertown • Maryland • 21620 • USA
800-344-3100 • 410-778-3100 (Outside U.S.A.) • Fax 410-778-6394
Visit us on the web at www.lamotte.com