

Certificate of Compliance for LaMotte Company


ISO 7027 Method

Testing performed on: 16, May 2005

ACMO Clear® co-polymer suspensions are tested for verification against formazin stock suspension from LaMotte Company. The AMCO Clear® calibration standard verification was performed in triplicate at five suspension levels against formazin per ISO 7027 section 6.2.4

Verification of proven AMCO Clear® polymer standards against a freshly prepared formazin concentrate should be performed once every 6 months to prove the measured average bias and precision of the AMCO Clear® polymer standards is no greater than the proven average bias and precision as determined by inter-laboratory studies (see annex A, ISO 7027)

Performed by: Michael McBride



Manager, AMCO Clear® Standards Division  
GFS Chemicals, Inc

Theoretical Values (FNU)	Formazin +/- 5% Lot # 244	AMCO Clear®		
		1	2	3
1.0	1.0	1.01	1.00	1.04
4.0	4.0	4.10	4.07	4.12
10	10	10.01	9.96	9.95
40	40	40.65	40.71	40.63
100	100	100	100	101